

tousimis Large Capacity Holders: 4", 3", 2" and 1"



Large Capacity Holder Compatibility

	# 8771A 1" Large Capacity Holder	# 8771B 2" Large Capacity Holder	# 8771C 3" Large Capacity Holder	# 8771D 4" Large Capacity Holder
SAMDRI®-PVT-3D	✓	X	X	X
SAMDRI®-795	✓	X	X	X
Autosamdri®-815	✓	X	X	X
Autosamdri®-815B	✓	✓	✓	✓
Autosamdri®-931 with 1.25" Chamber	✓	X	X	X
Autosamdri®-931 with 2.50" Chamber	✓	✓	X	X
Autosamdri®-931 with 3.40" Chamber	✓	✓	✓	X
Touchscreen Series C	✓	✓	✓	✓

Large Capacity Holders



4" Large Capacity Holder (Material: SS316) - Cat# 8771D
Interior 0.7" Ht x 4.0" I.D. (17.8 mm Ht x 101.6 mm I.D.)
Exterior 1.05" Ht x 4.4" O.D. (26.7 mm Ht x 111.8 mm O.D.)
Mesh Material: SS 316 / Mesh Size: 35x35 / Opening: 0.017" (0.43mm)



3" Large Capacity Holder (Material: SS316) - Cat# 8771C
Interior 0.7" Ht x 3.0" I.D. (17.8 mm Ht x 76.2 mm I.D.)
Exterior 1.05" Ht x 3.25" O.D. (26.7 mm Ht x 82.6 mm O.D.)
Mesh Material: SS 316 / Mesh Size: 35x35 / Opening: 0.017" (0.43mm)



2" Large Capacity Holder (Material: SS316) - Cat# 8771B
Interior 0.7" Ht x 2.0" I.D. (17.8 mm Ht x 50.8 mm I.D.)
Exterior 1.05" Ht x 2.25" O.D. (26.7 mm Ht x 57.2 mm O.D.)
Mesh Material: SS 316 / Mesh Size: 35x35 / Opening: 0.017" (0.43mm)



1" Large Capacity Holder (Material: SS316) - Cat# 8771A
Interior 0.7" Ht x 1.0" I.D. (17.8 mm Ht x 25.4 mm I.D.)
Exterior 1.05" Ht x 1.20" O.D. (26.7 mm Ht x 30.5 mm O.D.)
Mesh Material: SS 316 / Mesh Size: 35x35 / Opening: 0.017" (0.43mm)

4" Large Capacity Holder (Cat# 8771D)



3" Large Capacity Holder (Cat# 8771C)



2" Large Capacity Holder (Cat# 8771B)



1" Large Capacity Holder (Cat# 8771A)



Opening Large Capacity Holders



Large Capacity Holders

